

3D ex vivo imaging and characterization of vascular network in rodent heart



Tissues & Organs • Light Sheet • Advanced Media • Vascular Network • Cardiology

YOUR NEEDS

- Study of the impact of pathologies on the cardiac vascular network
- Preclinical study of treatment efficacy



General Procedure

Prior to sample collection by Imactiv-3D:

- In vivo labelling by infusion with a fluorescent lectin before euthanasia.
- Formalin fixation of extracted sample.

Image acquisition:


- Sample clearing.
- 3D light sheet fluorescence microscopy.
- Multi-position acquisition.

Image processing and analysis:

- Quantitative characterization of the vascular network:
 - Vessel segmentation.
 - Extraction of efficient volume.
 - Computation of parameters of interest: length and size of the vessels, density of the vascular network.
- 3D visualization with surface and volume rendering:
 - Reconstruction of the whole sample.
 - Advanced display using 3D animations.


OUR SOLUTIONS

- Light sheet microscopy and clearing to characterize the vascular network in 3D
- Automated 3D image processing for microvascular network quantification



Application example

in collaboration with SERVIER



- **Aim:** characterization of the microvascular network of a whole rat heart.
- **Technical specifications:**
 - Whole sample acquisition using 6 to 9 fields of view.
 - Raw data amounts to about 200 GB or 100 billion voxels.
 - Image processing takes 14 hours using a parallel architecture.
- Analysis of the microvascular network in 33 samples following this procedure exhibited an internal variation of only 3.5 % with respect to the mean.
- Assessment of micro-vascular alteration in pathological rat with comorbidities such as hypertension/diabetes (and effect of a potential treatment on vascular density).