

3D *ex vivo* imaging and characterization of neuromuscular junctions in rodent muscles

Tissues & Organs • Light Sheet • Advanced Media • Morphometric Quantification • Automatic Detection • Neurology



YOUR NEEDS

- Study of neuromuscular pathologies or atrophy models
- Preclinical evaluation of compound efficacy

OUR SOLUTIONS

- Light sheet microscopy and clearing to visualize neuromuscular junctions on whole muscle in 3D
- Automated image processing for quantification

General Procedure

Prior to sample collection by Imactiv-3D:


- *In vivo* labeling by infusion with a fluorescent toxin
- Formalin fixation of extracted sample

Image acquisition:


- Sample clearing
- 3D light sheet fluorescence microscopy
- Multi-position acquisition

Image processing and analysis:

- Quantitative characterization of all types of leg muscle:
 - Estimation of muscle volume
 - Absolute number and density of neuromuscular junctions
- 3D visualization with surface and volume rendering:
 - Reconstruction of the whole sample
 - Visualization of the spatial distribution of junctions
 - Advanced display using 3D animations



Application example: *in collaboration with* IPSEN Innovation for patient care



- Aim: comparison of a wide range of muscle types in the rat leg in terms of neuromuscular junction density
- Technical specifications:
 - Whole sample acquisition using 3 to 10 fields of view
 - Number of junctions ranges from 3 000 to 20 000 depending on the muscle type
 - Junction average size is about $60 \times 30 \times 10 \mu\text{m}^3$
- Analysis of the neuromuscular junction density following this procedure on 30 muscles produced interesting results among different types of muscle